TROLLEY MAINTENANCE LOG

rolley Serial #:						
Intake HEPA Change:/ Exhaust HEPA Changed:/						
		Air Exhaust				HEPA
		Pre-	& Hoses	Pre-	Trolly	Filters
Date	Initials	Filter	Cleaned		Certified	/Decon
enter date/initials and boxes						
		E2wks	E3mo	E12mo	E12mo	PRN

Trolley Air Exhaust Pre-Filter is changed every 2 weeks

- 1. With trolley running pre-filter is spritzed w/Sporicidin & bagged-out.
- 2.Trolley exhaust area is wiped w/Sporicidin & new filter installed.
- 3. Dirty exhaust filter bag is spritzed w/Sporicidin & removed for disposal

Trolley Rack Hoses are changed every 3 months

- 1. Hoses are removed & capped with glove do not create aerosols.
- 2.Hoses are spritzed w/Sporicidin & removed to dirty CW for cleaning/decon

Trolley Air Intake Pre-Filter is changed/sanitized every 12 months

Trolley Certification every 12 months

= 3 yrs/as needed - indicate intake (IT) or

exhaust (EX) HEPA